

Innovative, Intuitive, Flexible.

Luminex Flow Cytometry Solutions
with **Guava**[®] and **Amnis**[®] Systems

[Learn More >](#)



Luminex[®]
complexity simplified.

OMIP-028: Activation Panel for Rhesus Macaque NK Cell Subsets

Nicholas Pomplun,¹ Kim L. Weisgrau,¹ David T. Evans,^{1,2} Eva G. Rakasz^{1*}

¹Wisconsin National Primate Research Center, University Of Wisconsin-Madison, Wisconsin 53711

²Department of Pathology and Laboratory Medicine, University Of Wisconsin-Madison, Wisconsin 53711

Received 27 April 2015; Revised 3 July 2015; Accepted 10 July 2015

Grant sponsor: NIH, Grant numbers: 5P51OD011106-53; R01 AI095098

Additional Supporting Information may be found in the online version of this article.

*Correspondence to: Eva G. Rakasz, Ph.D., AIDS Vaccine Research Laboratory, 585 Science Drive, Madison, WI 53711, USA. E-mail: erakasz@primate.wisc.edu

Published online 28 July 2015 in Wiley Online Library (wileyonlinelibrary.com)

DOI: 10.1002/cyto.a.22727

© 2015 International Society for Advancement of Cytometry

• Key terms

NK cells; KIR3DL05; KIR3DL01; rhesus macaque

PURPOSE AND APPROPRIATE SAMPLE TYPES

THIS panel was developed to quantify natural killer (NK) cell subsets in Rhesus macaques (*Macaca mulatta*) during SIVmac239 infection induced pathogenesis. It includes markers to monitor changes in the activation/proliferation phenotype of up to 12 NK cell populations. The performance of the staining was tested on cryopreserved lymph node samples, and on fresh and cryopreserved peripheral blood mononuclear cells (PBMC) isolated from EDTA anti-coagulated blood. The panel can be used to characterize NK cells in a range of normal and pathologic conditions of this species and can be easily adapted to stain samples from various tissues.

BACKGROUND

NK cells are cytotoxic effector cells of the innate immune system. They are triggered by combined signals of activating and inhibitory receptors. In humans, NK cells are typically defined by the expression of CD56 (neural cell adhesion molecule—NCAM) and CD16 (Fc γ receptor IIIa) antigens with cells positive for both markers representing the overwhelming majority of the population. A subset of human NK cells express CD8 antigen. In the blood of healthy Rhesus macaques, CD16 and CD56 are co-expressed primarily by monocytes. The majority (>85%) of circulating NK cells are CD16 single positive. CD56 single positive, CD16/CD56 double positive, or double negative cells only represent a small fraction of the cells. In macaques, the total population of NK cells is defined by the expression of CD8 $\alpha\alpha$ chain and NKG2A/C antigens on CD3 negative lymphocytes (Fig. 1a) (1,2). The population can be further divided according to the expression of killer cell Ig-like (KIR) receptors. Currently, up to four NK cell populations can be delineated by two well-characterized KIRs—KIR3DL01 and KIR3DL05—in an individual animal (Fig. 1b) (3,4). Each subpopulation can be further divided into CD16+, CD56+, and CD16-CD56- subsets, altogether up to 12 subsets. Chronic simian immunodeficiency virus (SIV) infection induces a significant increase of the absolute number of NK cells and a major shift in the NK cell subset distribution (5). The change in the subset composition can be characterized by a pronounced increase of the CD16/CD56 double negative population. The CD56 single positive population always remains below 5%. The three subpopulations possess different Nkp46 (an NK cell stimulatory protein) expression profile (Fig. 1c), which shows animal-to-animal variation and may change depending on the stage of SIV-infection (5,6).

Activated/proliferating NK cells can express CD69, human leukocyte antigen-DR (HLA-DR), PD-1, and Ki-67 antigens to various extents (Fig. 1d). The combination of these markers defines 16 phenotypically distinct subsets probably representing

Table 1. Summary table for OMIP-028

Purpose	Characterize the activation and proliferation phenotype of various NK cell subsets
Species	Macaca mulatta
Cell types	Fresh PBMC and lymph node samples
Cross references	OMIP-007 (10)

different stages of maturation (Fig. 1e). With the present staining panel, we wanted to provide a tool to follow SIV infection-induced changes of NK cell subpopulations at a more refined level than previously described.

We had to work around two major limitations: scarcity of cross-reactive antibodies recognizing certain Rhesus macaque antigens, and limited number of reagents binding macaque KIR receptors. Mamu-KIR3DL01 is recognized only by the PE and APC conjugated version of monoclonal antibody NKVFS1. KIR3DL05 binds to a select number of Mamu-A1*00201 MHC-I allele restricted SIVmac239 epitope tetramers. Amongst these reagents the gag GY9 tetramer provided the best separation of the positive subsets (4). The GY9 tetramer is produced in BV421, PE, or APC conjugated forms by the NIH Tetramer Core Facility, all of which yield excellent staining. In Rhesus samples, the single available cross-reacting antibody—clone Z199—binds to both the activating NKG2C and the inhibitory NKG2A receptor (7). It can only be obtained in PE, PE-Cy7, and APC conjugated forms. Both APC and PE-Cy7 labeled versions provided excellent results. Of the two monoclonal antibody clones against the proliferation marker Ki-67 (B56 and Ki-67), only the Alexa Fluor 647 form of clone B56 yielded interpretable data. Finally, our only option of a cross-specific antibody for NKp46—clone BAB281—is available in PE, PE-Cy5, and PE-Cy7 conjugates. Our anchor reagents therefore

became KIR3DL05 (GY9) BV421, KIR3DL01 PE, NKp46 PE-Cy5, NKG2A/C PECy7, and Ki-67 Alexa Fluor 647.

There are multiple choices for cross-reacting antibody clones against the CD8 α chain. Since this protein is expressed abundantly on the cell surface, we selected a conjugate—the BV510 labeled SK1 clone—with good separation and least concern for broadening into other channels. In Rhesus macaques, the most frequently used clone against CD16 is the 3G8. It is available in 20 different fluorochrome-labeled versions from which we tested the V500, Pacific Blue, PE, PerCP-Cy5.5, and the BV711 varieties. All of the tested antibodies presented a clear separation between positive and negative populations. Since the Pacific Blue channel was reserved for KIR3DL05 BV421, the PE for KIR3DL01 PE and the PerCPCY5.5 channel for NKp46 PE-Cy5, we chose the BV711 conjugate. Tissue resident NK cell populations contain larger proportions of CD56 positive subsets than the circulating population (5). Whether this distribution shifts during acute or chronic infection remains to be determined. We tested three clones of CD56-specific antibodies conjugated either to FITC (clone B159), Vio-Bright FITC (clone AF12-7H3), PE-Cy7 (clone NCAM16.2), or Alexa Fluor 700 (clone B159). Although the PE-Cy7-labeled NCAM16.2 clone possessed the best stain index we chose the FITC labeled clone B159 for the final panel based on prioritization of detector blue A for NKG2A/C PE-Cy7 (Supporting Information Fig. 3).

Recently increased PD-1 expression on NK cells was implicated in better HIV-1 control (8). Therefore, we were interested to see whether similar correlation could be observed in SIV-infected Rhesus macaques. The sole cross-reactive PD-1-specific clone EH12.2H7 is available in more than ten fluorochrome-conjugated versions including six Brilliant Violet dye-conjugated forms. Since this antigen is not abundant on the cell surface it was an ideal candidate for the BV605 conjugate. This fluorochrome has a second emission peak at 650 so there is a

Table 2. Reagents used in OMIP-028

SPECIFICITY	CLONE	FLUOROCROME	PURPOSE
Live/Dead	N/A	Near infrared	Viability
CD45	D058-1283	BV786 ^a	hematopoietic cell lineage
CD3	SP34-2	Alexa Fluor 700	Exclusion
CD20	2H7	Alexa Fluor 700	
CD14	M5E2	Alexa Fluor 700	
CD8	SK1	BV510	NK subsets
CD16	3G8	BV711	
CD56	B159	FITC	
NKG2A/C	Z199	PE-Cy7	
KIR3DL01	NKVFS1	PE	
KIR3DL05	GY9 tetramer	BV421	
NKp46	BAB281	PE-Cy5	NK cell differentiation
CD279 (PD-1)	EH12.2H7	BV605	
CD69	TP1.55.3	ECD ^b (PE-Texas Red)	Activation
HLA-DR	L243	BV650	
Ki-67	B56	Alexa Fluor 647	Proliferation

^aBV = Brilliant Violet.

^bECD = energy coupled dye (PE-Texas Red).

possibility of significant spillover into the Violet C channel. The low intensity staining of the PD-1 BV605 conjugate provided us with an acceptable staining pattern with no fluorescence spillover. There are four cross-reactive monoclonal antibodies against the very early activation antigen CD69. We tested the PE-CF594 labeled FN50 and the ECD labeled TP1.55.3 clones. Both of them yielded similar results regarding the stain index and frequency of positive events. Finally, HLA-DR was shown to be present on a subset of NK cells characterized by enhanced degranulation (9). This antigen is expressed approximately one log lower on NK cells compared to B cells or monocytes. There is an abundance of cross-reactive HLA-DR-specific monoclonal antibodies in more than 50 fluorochrome-conjugated forms. We evaluated clone L243 in APC-Cy7 and BV650 labeled form. Both versions showed clear separation of HLA-DR positive NK cell populations. Since we reserved detector Red A for the live/dead discrimination we opted for the BV650 conjugate.

ACKNOWLEDGMENTS

Animals were handled in accordance with the standards of the American Association for the Accreditation of Laboratory Animal Care (AAALAC).

LITERATURE CITED

1. Hermes M, Albrecht C, Schrod A, Brameier M, Walter L. Expression patterns of killer cell immunoglobulin-like receptors (KIR) of NK-cell and T-cell subsets in Old World monkeys. *PLoS One* 2013;8:e64936.
2. Biassoni R, Fogli M, Cantoni C, Costa P, Conte R, Koopman G, Cafaro A, Ensoli B, Moretta A, Moretta L, et al. Molecular and functional characterization of NKG2D, NKp80, and NKG2C triggering NK cell receptors in rhesus and cynomolgus macaques: Monitoring of NK cell function during simian HIV infection. *J Immunol* 2005; 174:5695–5705.
3. Colantonio AD, Bimber BN, Neidermyer WJ Jr, Reeves RK, Alter G, Altfeld M, Johnson RP, Carrington M, O'Connor DH, Evans DT. KIR polymorphisms modulate peptide-dependent binding to an MHC class I ligand with a Bw6 motif. *PLoS Pathog* 2011;7:e1001316.
4. Schafer JL, Colantonio AD, Neidermyer WJ, Dudley DM, Connole M, O'Connor DH, Evans DT. KIR3DL01 recognition of Bw4 ligands in the rhesus macaque: Maintenance of Bw4 specificity since the divergence of apes and Old World monkeys. *J Immunol* 2014;192:1907–1917.
5. Reeves RK, Gillis J, Wong FE, Yu Y, Connole M, Johnson RP. CD16- natural killer cells: Enrichment in mucosal and secondary lymphoid tissues and altered function during chronic SIV infection. *Blood* 2010;115:4439–4446.
6. Mandelboim O, Lieberman N, Lev M, Paul L, Arnon TI, Bushkin Y, Davis DM, Strominger JL, Yewdell JW, Porgador A. Recognition of haemagglutinins on virus-infected cells by NKp46 activates lysis by human NK cells. *Nature* 2001;409:1055–1060.
7. LaBonte ML, Choi EI, Letvin NL. Molecular determinants regulating the pairing of NKG2 molecules with CD94 for cell surface heterodimer expression. *J Immunol* 2004;172:6902–6912.
8. Taborada N, Hernandez JC, Laiola J, Juno JA, Kimani J, Rugeles MT, Fowke KR. Low expression of activation and inhibitory molecules on NK cells and CD4+ T cells is associated with viral control. *AIDS Res Hum Retroviruses* 2015;31:636–640.
9. Evans JH, Horowitz A, Mehrabi M, Wise EL, Pease JE, Riley EM, Davis DM. A distinct subset of human NK cells expressing HLA-DR expand in response to IL-2 and can aid immune responses to BCG. *Eur J Immunol* 2011;41:1924–1933.
10. Currier JR, Eller MA. OMIP-007: Phenotypic analysis of human natural killer cells. *Cytometry A* 2012;81A:447–449.